



This manual covers the use and maintenance for the **CS100** controller installed on **the LS SERIES Cryogenic Storage System**. Please read all documentation before using this equipment and all personnel who use this equipment should be trained. In addition, all service and maintenance on this equipment should be performed by an authorized distributor, service company or technician.

Failure to follow instructions in this manual may result in damage to equipment, poor or unexpected performance or injury to personnel operating the equipment.

Material in this manual is for information purposes only. The contents and the product it describes are subject to change without notice. ICBiomedical, LLC, makes no representations or warranties with respect to this manual. In no event shall ICBiomedical, LLC, be held liable for any damages, direct or incidental arising out of or related to the use of this manual.





CAUTION - SAFETY FIRST!

- **REVIEW AND UNDERSTAND ALL SAFETY PROCEDURES IN THE "HANDLE WITH CARE" BOOKLET FORM # TW-10 P/N 7950-8052 BEFORE ATTEMPTING TO INSTALL, OPERATE OR PERFORM MAINTENANCE ON THIS LN2 CONTROLLER**
- **DO NOT ATTEMPT TO USE OR MAINTAIN ANY LIQUID NITROGEN FREEZER UNTIL YOU READ AND UNDERSTAND THESE INSTRUCTIONS**
- **DO NOT PERMIT UNTRAINED PERSONS TO USE OR MAINTAIN THIS UNIT**
- **IF YOU DO NOT FULLY UNDERSTAND THESE INSTRUCTIONS, CONTACT YOUR SUPPLIER FOR FURTHER INFORMATION**
- **BEFORE ATTEMPTING TO OPERATE THIS CONTROLLER WITH ANY WORTHINGTON INDUSTRIES™ LS SERIES UNITS, YOU MUST READ THE SEPARATE OPERATING AND SAFETY MANUAL PROVIDED WITH THAT WORTHINGTON INDUSTRIES UNIT**

Table of Contents

Symbols.....	5
Safety	5
Electromagnetic Compatibility (EMC)	6
Operation	6
Operating Environment.....	6
General Equipment Description	7
Main Control Front Panel	7
Main Control Back Panel	8
Features	8
LN2 Level Measurement and Control	8
Temperature Monitoring.....	9
Alarm Conditions	9
Operation Data	10
Communications	10
Control Setting Adjustments	11
Level.....	11
Menu System.....	11
Thermocouple Enable	11
Set Date/Time	12
Low Supply Alarm	12
High Temp Alarm	12
Temp Alarm Delay	12
Audible Alarm Retrigger	12
Remote Alarm Delay	12
Display Brightness.....	13
About	13
Menu System.....	13
Factory Defaults.....	13
Validation	13
Maintenance	14
Troubleshooting	14
Product Information	14
Specifications.....	14
Interconnection Block Diagram.....	16
Replacement Parts.....	17
Appendix	18
Certifications & Listings.....	18
EN Compliance Tables	18
Declaration of Conformity	22
Warranty	22
Installation & Setup	23
Temperature Calibration	25
External Connector Ratings	26
Labels	26

Symbols

WARNING! Avoid injury	
WARNING! Low temperature	
Read owner's manual before use.	
Waste electrical and electronic equipment directive	

Safety



Please use proper handling procedures when working with Liquid Nitrogen (LN₂). It is an extremely cold refrigerant and can cause serious injury if not handled properly. Also, Nitrogen gas can cause asphyxiation in a confined area so be sure that all confined areas have adequate ventilation. For detailed information on the handling of cryogenic liquids, refer to the publication: P12 "Safe Handling of Cryogenic Liquids" available from the Compressed Gas Association Inc., 1235 Jefferson Davis Highway, Arlington, VA 22202.
www.cganet.com

- This product is not intended for a life support function
- This product is intended to be used in offices/ clinics as well as hospitals
- This product has no Radio Transmitter (Intentional Radiator) functions
- This product is not intended for electromagnetic shielded rooms only
- This product does not intentionally apply RF energy for its function
- This product does not intentionally receive RF energy for its function
- This product is not a large permanently- installed product that is part of a system

Sample safety is of vital importance when using a Cryo-Storage system. Many times, the product stored represents many years of research and development or the samples may represent someone's hopes and dreams to have children or to recover from an illness. It is of critical importance to avoid problems through due diligence. The control system is a single aspect of the Cryo-Storage system required to ensure proper sample safety.

Additional aspects include:

- **Planning:** Plans should be implemented such that no confusion exists regarding the stored product, the equipment or the personnel.
- **Procedures:** Procedures should be implemented to best protect the samples. This would include proper handling so that sample viability is never compromised. This should include a manual dipstick of the LN2 level on a daily basis.
- **Training:** Personnel should be trained on a periodic basis including refresher courses and critical evaluation of planning and procedures.
- **Documentation:** All aspects of the system should be thoroughly documented with special emphasis given to deviations from normal operations.
- **Maintenance:** Equipment should be maintained to top condition through a preventive maintenance program.
- **Redundancy:** Redundant systems should be used to ensure that a single point of failure cannot compromise sample safety. This would include additional sensors to monitor critical parameters in the Cryo-Storage system. It would also include separating samples and storing in multiple Cryo-Storage vessels.
- **Risk Management:** Risks should be evaluated and minimized to ensure the safety of the samples.

Responsibility for sample safety ultimately lies with the owners and users of the Cryo-Storage equipment.

Electromagnetic Compatibility (EMC)

Although this equipment conforms to the intent of the 2004/108/EC EMC Directive, all medical equipment may produce electromagnetic interference or be susceptible to electromagnetic interference. The following are guidance and manufacturer's declarations regarding EMC for the Model CS100.

The Model CS100 needs special precautions regarding EMC and needs to be installed and put into service according to the RMA information provided in the following pages.

As with all electrical medical equipment, this equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures such as re-orienting or relocating the Model CS100 FOR LS SERIES unit or shielding the location.

WARNING:
**This equipment is intended
for use by healthcare
professionals.**

Portable and Mobile RF communications equipment can affect the performance of the Model CS100. Please use the guidelines and recommendations specified in Tables 4 and 6.

Other Medical Equipment or Systems can produce electromagnetic emissions and therefore can interfere with the functionality of the Model CS100. Care should not be used when operating the Model CS100 adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Model CS100 should initially be observed to verify normal operation in the configuration in which it will be used.

The electrical cables, external power supplies and accessories listed or referenced in this manual have been shown to comply with the test requirements listed in the following tables. Care should be taken to use only manufacturer-recommended cables, power supplies and electrical accessories with the Model CS100. If a third-party supplier offers cables, external power supplies and electrical accessories for use with the Model CS100 and they are not listed or referenced in this manual, it is the responsibility of that third-party supplier to determine compliance with the standards and tests in the following tables.

The use of electrical cables and accessories other than those specified in this manual or referenced documents may result in increased electromagnetic emissions from the Model CS100 or decreased electromagnetic immunity of the Model CS100.

Operation

Operating Environment

The CS100 is designed to be operated at normal room temperatures (60°F to 80°F, 15°C to 27°C) and a relative humidity level near 45%. The humidity level should be maintained such that the electronics are not exposed to condensation.

The LS SERIES Cryostorage system should be positioned such that all sides of the unit are easily accessible and the user can easily connect/disconnect the power cord from the wall socket.

Operator Skills, Training, and Knowledge

The personnel using the CS100 for LS Series and the Cryo-Storage freezer should have a basic understanding of the safety procedures involved when working with LN₂. In addition, the operating personnel should understand the messages and error conditions associated with the control system and should request support if questions arise.

Training and support are available through the distributor who sold the equipment as well as through the manufacturer of the equipment.

CE Control System Components

The CS100 for LS Series for the Cryo-Storage system consists of the following components.


- Control/Display Module
- Wiring Harness Assembly
- Power Supply
- Thermocouple Assembly
- Sensor Assembly
- Cryogenic Solenoid Valve
- Remote Alarm Plug
- Freezeguard Sensor

Some configurations may differ depending on the Cryo-Storage system being used.

The CS100 is assembled onto the freezer at the factory and completely tested prior to re-boxing for shipment. Refer to the Quick Start Guide for freezer setup. Start operation of the control system by plugging the power cord into the wall outlet. The CS100 will go through a short startup routine and then start operation. For information relating to assembly of the control components and connection information, please refer to Appendix A.


General Equipment Description

Main Control Front Panel



Operation

Main Control Under Panel



- A. Local Cryowire port
- B. Lid Switch
- C. Remote Alarm
- D. Solenoid Valve and Freezeguard Sensor
- E. Power (12VDC)
- F. Level Sensor
- G. Temperature Sensor

Features

The CS100 Automatically maintains the Liquid Nitrogen (LN₂) level and monitors temperature in the Cryogenic freezer. In addition, operational conditions are monitored, and an alarm is triggered if necessary. Audit and operations data is stored in memory on the control board.

LN₂ Level Measurement and Control

The LN₂ level gauge on the face of the control indicated level by lighting the appropriate LED if a thermistor is submerged in LN₂.

The CS100 uses thermistors to measure the LN₂ level within the vessel. A thermistor is a thermal resistor and its resistance changes as the temperature changes. When a thermistor is submerged in LN₂, its resistance will be significantly greater than its resistance at room temperature. The control can detect this resistance change and determine the level of the LN₂ within the freezer. The CS100 uses 4 thermistors to maintain the level. These thermistors correspond to Low (Low Alarm), Normal (Start Fill), Normal (Stop Fill) and High (High Alarm). When the LN₂ level drops below the Start Fill thermistor, the control opens a solenoid valve allowing LN₂ to enter the vessel. This continues until the Stop Fill thermistor is submerged in LN₂ at which point the solenoid valve is closed preventing the flow of additional LN₂ into the vessel. The High Alarm thermistor is located 1" above the Stop Fill thermistor and provides a safety. Likewise, the Low Alarm thermistor is located 1" below the Start Fill thermistor and provides a safety. Factory settings provide a 2" range between the start fill and stop fill thermistors.





The LN₂ level gauge on the face of the control indicates level by lighting the appropriate LED if a thermistor is submerged in LN₂. In the picture, the Low Alarm thermistor, the Start Fill thermistor and the Stop Fill thermistor are all submerged in LN₂.

Temperature Monitoring

The CS100 uses a Type T Thermocouple to determine the temperature in the vessel. The thermocouple is installed near the top of the LS SERIES freezer and the temperatures are measured at that point. Multiple temperature values are taken per second and then averaged.

The control provides a High Temperature Alarm which can be adjusted by the user. If the temperature exceeds the temperature alarm set point, the display flashes and an audible alarm is triggered.

Please see Appendix B for instructions on calibrating the temperature.



Alarm Conditions

The CS100 monitors a number of conditions and provides an alarm if a problem is detected. The alarms are listed below:

Low Level Alarm

LN₂ level is too low. Thermistor #1 on the sensor assembly is not submerged in LN₂.

High Level Alarm

LN₂ level is too high. Thermistor #4 on the sensor assembly is submerged in LN₂.

Sensor Fault Alarm

A problem exists with the level sensor. The control detects an open sensor circuit, meaning that the sensor is unplugged or the sensor assembly has been damaged.

High Temp Alarm

The temperature detected exceeds the high temperature alarm setting.

Thermocouple Open Alarm

A problem exists with the temperature sensor (thermocouple). The control detects an open circuit, meaning that the sensor is unplugged or the sensor assembly has been damaged.

Power Failure

No Power.

Low LN2 Supply Alarm

A problem may exist with the LN2 supply connected to the freezer. This alarm occurs if the freezer does not fill within the designated amount of time determined by the setting on the control. This may occur for a number of reasons including an empty supply cylinder, low pressure in the supply cylinder or a closed shut off valve.

Temp Alarm Delay

This is the amount of time after a warm temperature is detected before the control goes into alarm.

Audible Alarm Retrigger

The audible alarm is retriggered if the error condition that caused it is not corrected. The retrigger time can be adjusted by the user.

Remote Alarm Delay

The control provides a relay to provide an external signal that an alarm condition has occurred. The user can set the remote alarm timer that determines the amount of time an error must be active before the relay is triggered.

All alarms include the following:

- The LED array screen flashes to signal an error condition
- An audible tone sounds
- The error detected is displayed and scrolled on the screen
- The remote alarm relay changes state to provide a dry contact output signal

Operation Data

The CS100 stores data related to the operation of the Cryostorage freezer. This data includes date, time, LN2 level, temperature, system events and error conditions. This data can be useful for audit purposes, operation analysis and preventive maintenance.

Communications

The CS100 has been designed with advanced communications capabilities. This allows for the transfer of data out of the control where the data can more easily be used. Please check with your supplier for available protocols and compatible products.

Control Setting Adjustments











Level

The standard sensor assembly that is installed on a freezer consists of a circuit board with thermistors installed. The assembly has a fixed range of 4". The level can be changed by either raising or lowering the circuit board within the freezer.

If a wider range is needed between the start fill and stop fill thermistor, please contact your supplier

Menu System

Some control settings can be changed through the menu system.

Enter the menu system by pressing:	 
Move down through the menu system by pressing:	 
Move up through the menu system by pressing:	 
Select a menu choice or lock in a setting by pressing:	 
Back out of the menu system by pressing:	 

When changing settings, single button presses will increment/decrement a value one step at a time. Pressing and holding a button will allow for rapid change of a setting.




The control will continue to monitor all sensors and conditions while the user accesses the menu system. If no activity is detected for 30 seconds, the control will automatically revert to the main operational screen.

Thermocouple Enable

The temperature sensor (thermocouple) can be turned on / off through the menu system. If the thermocouple is turned on, the temperature is displayed on the LED display and a high temperature alarm will occur if the temperature exceeds the user-defined high temperature alarm setting. If the temperature sensor is turned off, no temperature is displayed and no high temperature alarm can occur.

Set Date/Time

The date and time can be set through the menu system. The date and time is used to provide a date stamp for data that is collected within the control. Use the appropriate buttons to adjust the settings and lock in the values.

Select Date/Time and lock in values:	
Increment Values:	
Decrement values:	

Low Supply Alarm

The Low Supply Alarm timer can be changed through the menu system and can range from 30 to 120 minutes. This alarm occurs if the filling operation takes too long to complete. If this alarm occurs the supply should be checked to insure that the supply valve on turned on, a sufficient supply of LN2 is available in the supply source and the supply source pressure is adequate.

High Temp Alarm

The High Temp Alarm setting determines the temperature at which the control will indicate a high temperature alarm. This is accessed through the menu system and values can range from 0°C to -196°C.

Temp Alarm Delay

The Temp Alarm Delay setting determines the length of time after the control detects a warm temperature before it goes into alarm. This is accessed through the menu system and values can range from 0 to 60 minutes.

Audible Alarm Retrigger

The Audible Alarm Retrigger setting determines the length of time before the audible alarm reoccurs after it has been muted. This only occurs if the error condition has not been corrected. This is accessed through the menu system and values range from 0 to 720 minutes in 10 minutes increments.

Remote Alarm Delay

The Remote Alarm Delay setting determines the length of time before the remote alarm relay is triggered after an error condition occurs. This is accessed through the menu system and values range from 0 to 720 minutes in 10 minute increments.

Display Brightness

The Display Brightness is the setting that determines the lamination intensity of the display. This is accessed through the menu system and values range from 1 to 15.

About

This provides information about the control including the firmware version, the control serial number, and contact information. This is accessed through the menu system.

Factory Defaults

- Thermocouple..... On
- Date/Time Factory Set for Easter Time
- Low Supply Alarm 30 Minutes
- High Temp Alarm -100°C
- Temp Alarm Delay 0 Minutes
- Audible Repeat Delay 30 Minutes
- Remote Alarm Delay..... 30 Minutes
- Auto Defog..... 30 Seconds
- Lid Open Alarm 600 Seconds (10 minutes)
- Lid Switch..... Off
- Quick Chill 30 Seconds
- Display Brightness 7

Validation

Some organizations require that equipment be validated periodically. If information is needed on the proper techniques to validate this equipment, please contact your supplier.

Maintenance

To ensure proper operation and maintain excellent performance of the LS SERIES, a regular maintenance schedule should be followed for the CS100. This would include the following:

CS100 Control System	Examine for exposure to moisture, wear and tear, connector problems, and damage to the faceplate or buttons. In addition, periodic firmware updates may be important.
Solenoid Valve	Examine wires and connector for damage. Replace every 3 years.
Level Sensor Assembly	Examine wires and connector for damage. Replace every 3 years.
Thermocouple Assembly	Examine wires and connector for damage.
Power Supply	Examine for damage to power supply and power cords.

In addition, inspection and preventive maintenance should also be performed on the freezer and its mechanical parts. Refer to owner's manual for details.

If any intermittent operation with the Cryostorage System is observed or suspected, it should be investigated and remedied immediately even if this falls outside of the normal maintenance schedule.

Troubleshooting

If the LS SERIES freezer experiences problems or appears to not be operating at optimum efficiency, please contact your supplier for assistance. The CS100 has incorporated state of the art diagnostic tools to assist in the identification and correction of any issues that may arise.


Product Information

Specifications	
Control Type	LN ₂ Level Control
Level Measurement	
Sensor Type	4-Thermistor Fixed 4- Thermistor Adjustable
Range	Low, Normal, High
Redundancy	Multiple discrete points
Temperature Measurement Sensor Type	Type T Thermocouple
Accuracy	1° or 1.5% of reading, whichever is greater
Resolution	.1°C
Number of Channels	1
Temperature Display Units	°C
Electrical	
Input Voltage	100-240 VAC
Input Current (max)	1.75 A
Input Current (continuous)	.5 A
Power Consumption (max)	21 W
Power Consumption (continuous)	6 W
Frequency	50/60 Hz
Output	12 VDC
Control Input Voltage	12 VDC
Power cord	Available for all countries
Solenoid Valve	
Input Voltage	12 VDC
Input Current	.96 amps
Communications	
Protocol	CAN
Number of Communication Ports	1
User Interface	
Display Type	LED Array
Buttons	4
Level Display	4 LED's
Filling	1 LED

Specifications Continued	
Control Tests	
Power Up Self-Test	Control system check
Thermistor Status	Yes
Alarms	
Low Level Alarm	Always Enabled
High Level Alarm	Always Enabled
Sensor Error Alarm	Always Enabled
High Temperature	Alarm Programmable
Thermocouple Open Alarm	Always Enabled
Remote Alarm	Always Enabled
Power Failure (Remote Only)	Always Enabled
Low Supply	Alarm Programmable
Audible Alarm	Always Enabled (Re-trigger Programmable)
Visual Alarm Indicator	Always Enabled
Buttons	
Fill (up arrow)	Open Valve (Menu up)
Stop Fill (down arrow)	Close Valve (Menu down)
Menu (Enter)	Access Menu (Save setting)
Mute (Back)	Silence Audible (Exit menu level)
Data Collection	
Temperature	Yes
Level	Yes
Alarms	Yes
Memory	4Mb
Dimensions	
Width	9.5" (241 mm)
Height	2.0" (51 mm)
Depth	1.31" (33 mm)
Weight	0.625 lbs. (0.28 kg)


Operation

Interconnection Block Diagram



Parts List

The parts and components listed below have been specified and tested for use with the CS100. These are not user serviceable parts. Replacement parts should be obtained from your distributor or supplier.



BILL OF MATERIALS			
ITEM NO.	PART OR CODE NO.	QTY	DESCRIPTION
1	5140-1310	1	CONTROLLER, LS6000-CS100
2	RS60-9C22	2	BRACKET, CONTROLLER
3	6999-9041	1	12 VDC CRYOGENIC SOLENOID VALVE CS200
4	6816-0025	1	STREET TEE, 1/4", BRASS
5	6719-9037	1	REDUCING HEX NIPPLE, 3/8 NPT X 1/4 NPT, BRASS
6	7631-1080	1	STRAINER, 3/8 NPT, BRASS BODY, SS ELEMENT
7	7355-4712	1	CONNECTOR CGA-295, 1/2 ODT-45 FLARE X 3/8 MNPT, BRASS
8	45701960	1	CONNECTOR, MALE, 3/8" ODT X 1/4" MNPT, BRASS
9	RS60-9C37	1	RELIEF VALVE TUBE
10	45701970	1	FEMALE CONNECTOR, .38 ODT COMPRESSION X .25 FPT, BRASS
11	6913-9077	1	RELIEF VALVE, 1/4" MNPT, 100 PSIG, BRASS
12	45251806	1	ELBOW, 90 DEG., 1/4" MNPT X 3/8" ODT, S.S.
13	LS60-9C33	1	SENSOR & FILL TUBE ASSEMBLY, LS6000
14	6460-0884	8	WASHER, FLAT, 1/4", STN STL
15	41431200	4	CAP SCREW, HEX HEAD, .25-20 UNC X 1.0 LG, SS
16	6460-2025	4	WASHER, LOCK, 1/4", STN STL
17	38221798	4	ELASTIC STOP NUT, .25-20 UNC, BRASS
18	6460-7090	2	WASHER, EXTERNAL TOOTH, #10
19	6160-4905	2	SCREW, MACHINE, HEX HEAD, 10-32 UNC X 0.44 LG
20	9365-0001	2 FT	CABLE WRAP, SPIRAL CUT, TEFLON, 0.25
21	9365-0002	4 FT	CABLE WRAP, SPIRAL CUT, TEFLON, 0.50
22	7601-1025	4	FASTENER, SNAP-ON, M-240
23	7200-4010	1	BAG, POLY, 0.0015 X 3 X 4
24	7200-4016	1	BAG, POLY, 0.004 X 8 X 12
25	7200-4055	1	BAG, POLY, 0.004 X 32 X 42
26	7200-4043	1	BAG, POLY, 0.004 X 6 X 9
27	6033100	1	INSTALLATION INSTRUCTIONS FOR LS6000 WITH CS100 CONTROLLER
28	6460-3005	2	WASHER, FLAT, #6, STAINLESS STEEL
29	3701-9106	1	CARTON, SHIPPING, 8" WIDE X 8" TALL X 40" LONG
30	RS60-8C30	1	TRANSFORMER BRACKET, LS6000-CS100
31	LA80-9C51	1	TRANSFORMER COVER EL-V LABS CS
32	5140-1229	1	12 VDC TRANSFORMER, NO PLUG
33	6160-9920	4	SCREW, MACHINE, TRUSS HEAD, 10-32 UNC X 1/2" LG
34	6160-3900	2	SCREW, ROUND HEAD, 8-32 UNC X 1/2 LG, SS
35	5140-1281	1	4 THERMISTOR SENSOR ASSEMBLY, CS100
36	5140-1289	1	SOLENOID VALVE HARNESS ASSEMBLY, LS6000-CS100
37	5140-1219	1	THERMOCOUPLE
38	373491	1	MANUAL, LS-SERIES CS100 CTRL
39	5140-1230	1	POWER CORD, 100-240 VAC, U S A
40	7200-4312	2	BAG, POLY, 0.004 X 12 X 14

NOTE:

- FUNCTION TEST PER FORM NO 330.
- USE TEFLON TAPE ON ALL THREADED JOINTS.
- PLACE ALL FASTENERS IN BAG (ITEM 23) AND HEAT SEAL CLOSED. REPACK CONTROLLER IN ITS BOX. AFTER ASSEMBLING ITEMS 9, 10, AND 11, PLACE IN BAG (ITEM 24) WITH CONTROLLER BRACKETS (ITEM 2), AND TRANSFORMER BRACKET (30) THEN PLACE BAG OF FASTENERS IN AND HEAT SEAL CLOSED.
- INSTALL SENSOR (ITEM 35) AND THERMOCOUPLE (ITEM 37) INTO FILL/SENSOR TUBE ASSEMBLY USING CABLE WRAP (ITEM 21). PLACE ITEM 20 IN BAG (ITEM 26) AND HEAT SEAL CLOSED. (REFER TO LS6000-CS100 DRAWING FOR ASSEMBLY INSTRUCTIONS OF SENSOR TUBE TO FILL TUBE).
- AFTER ASSEMBLING PLUMBING COMPONENTS (ITEMS 3, 4, 5, 6, 7, 8, 12, 31 AND 32), PLACE IN BAG (ITEM 27) AND HEAT SEAL CLOSED.
- PLACE COPY OF INSTALLATION INSTRUCTIONS (ITEM 27), AND MANUAL (ITEM 38) IN BAG (ITEM 40).
- PLACE POWER CORD (ITEM 39) IN BAG (ITEM 40).
- PLACE ALL BAGGED ITEMS IN BOX (ITEM 29) AND FILL WITH PACKING PEANUTS. PLACE IN BOX AND SEAL BOX CLOSED WITH FILAMENT SHIPPING TAPE. PLACE PART NUMBER LABEL ON ALL FOUR SIDES OF THE BOX.

Certifications & Listings

This product complies with the following standards and directives:

- IEC 61010-1
- UL 61010-1
- CAN/CSAC22.2#61010-1
- RoHS Directive
- ETL/cETL Listing for North America
- CE marked to the Low Voltage Directive

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The CS100 is intended for use in the electromagnetic environment specified below. The customer or the end user of the CS100 should assume that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions - CISPR 11 (Radiated & Conducted)	Group 1	The CS100 FOR LS SERIES uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions - CISPR 11 (Radiated & Conducted)	Class A	The CS100 FOR LS SERIES is suitable for use in all commercial establishments other than domestic, and many be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded. Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the CS100 FOR LS SERIES or shielding the location.
Harmonic Emissions EN/IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker Emissions EN/IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

Recommended separation distances between portable and mobile RF communications equipment and the CS100

The CS100 FOR LS SERIES is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CS100 FOR LS SERIES can help prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communications equipment (transmitters) and the CS100 FOR LS SERIES as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter in watts (W)	Separation distance according to frequency of transmitter in meters (m)		
	150kHz to 80MHz d = 1.2	80MHz to 800 MHz d = 1.2	800MHz to 2.5GHz d = 2.3
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1.0	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

EC Declaration of Conformity

The undersigned representing the manufacturer:

Pacer Digital Systems, Inc.

Attn: Kevin Oeff

8658 Castle Park Drive, Suite 103

Indianapolis, IN 46256, USA

Herewith declared that the Product: LN₂ Level Control for Cryostorage System

Model/Type ref.: CS CONTROL SYSTEM

is in conformity with the essential requirements of the following EC Directives when subject to the correct installation, maintenance and use conforming to the intended purpose of the product. The conformance in regard to the applicable regulations and standards is consistent with operations and maintenance in the Owner's manual. This declaration has been issued under the sole responsibility of the manufacturer. The object of the declaration is in conformity with the relevant Union harmonization Legislation.

2014/35/EU Low Voltage Directive

2014/30/EU EMC Directive

and that the Standards and/or technical specifications referenced below have been applied:

- CENELEC EN 60601-1:2006 (Third Edition): Medical Electrical Equipment – General Requirements for basic safety and essential performance.
- UL 61010-1:2012 Ed.3+R:29Apr2016]: Safety Requirements For Electrical Equipment For Measurement, Control, And Laboratory Use.
- CSA C22.2#61010-1-12:2012 Ed.3+U1;U2]: Safety Requirements For Electrical Equipment For Measurement, Control, And Laboratory Use.
- IEC 60601-1-2: 2007 (Third Edition) Medical Electrical Equipment – General Requirements for basic safety and essential performance – Collateral standard : Electromagnetic Compatibility.
- IEC/CISPR 11:2009+A1:2010 Radiated & Conducted Emissions.
- IEC 61000-3-2:2005+A1:2008+A2:2009. Harmonics.
- IEC 61000-3-3:2008. Flicker

Year of CE Marking: 2012

The Technical Construction File required by this Directive is maintained at the corporate headquarters of Pacer Digital Systems, Inc., 8658 Castle Park Drive, Suite 103, Indianapolis, Indiana, USA.

Manufacturer: Pacer Digital Systems, Inc.

Signature:

Position: President

Date: 23 Dec 2021

Place: Indianapolis, IN USA

Warranty

ICBiomedical, LLC, warrants that each of its electronic control products will be free from defects in material and workmanship, in the normal service for which the product was manufactured, for a period of two years from the date of purchase. Parts and components used in conjunction with electronic controllers and displays are warranted to be free from defects in material and workmanship, in the normal service for which the product was manufactured, for a period of ninety (90) days from the date of purchase. ICBiomedical, LLC, at its option discretion, will either repair or replace any item covered under this warranty.

This warranty is void if the product is used for any other purpose than that for which it was designed, including but not limited to connection with third party systems. This warranty is also void if the product is in any way altered or repaired by others. ICBiomedical, LLC, shall not be liable under this warranty, or otherwise, for defects caused by negligence, abuse or misuse of this product, corrosion, fire or the effects of normal wear.

The remedies set forth herein are exclusive. ICBiomedical, LLC, shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products, resulting from the delivery, use or failure of the product or for any other cause. By accepting delivery of the product, the purchaser acknowledges that this limitation of remedies is reasonable and enforceable. In no case shall ICBiomedical, LLC, liability exceed the purchase price for the product.

Installation & Setup

The CS100 consists of the following components.

- Main Control / Display
- Power Supply
- Thermocouple Assembly
- Sensor Assembly
- Cryogenic Solenoid Valve
- Remote Alarm Plug

Connect the wiring harness assembly to the main control/display. The connector is keyed and can only be plugged in one way.

Connect the thermocouple plug of the harness assembly into the control at the thermocouple plug. This is located on the right edge of the control if looking at it from the front. One of the blades on the plug is slightly wider insuring that it is plugged in correctly. Copper blade should plug into copper colored plug. On the other end of the harness, connect the Thermocouple assembly into the plug insuring that copper colored blade is plugged into copper colored contact on receptacle.

Connect the level sensor assembly into the panel at the end of the level sensor harness assembly. The connector is a round locking connector that is keyed so it can only be plugged in one way.

Connect the solenoid valve into the panel at the end of the wiring harness assembly. The connector is a 4 pin locking connector that is keyed so it can only be plugged in one way.

Connect the Remote Alarm plug into the panel at the end of the wiring harness assembly. This is keyed so that it can only be plugged in one way.

Finally, connect the barrel plug of the power supply into the receptacle on the back of the controller panel.

Temperature Calibration

The CS100 control system uses a type T thermocouple to measure temperature within the Cryostorage vessel. The temperature curve for a thermocouple is nonlinear so it is important that the CS100 have a good calibration to provide accurate temperature readings. The temperature is traceable to the National Institute of Standards and Technology (NIST) ITS-90 Thermocouple Database.

There are three important reference points needed for calibration:

- Ambient Temperature..... **
- Ice water 0°C
- Liquid Nitrogen (LN₂) -196°C


**The ambient temperature is the temperature measured inside the control box and is used to provide temperature compensation adjustment. This is calibrated at the factory and should not be adjusted in the field.

Ice water and LN₂ provide the reference points on the temperature curve. If these two points are calibrated correctly then all other points (temperatures) on the curve are correct. The accuracy of a thermocouple is +/- 1°C or +/- 1.5% of the reading, whichever is greater.

To check a calibration, dip the thermocouple in ice water and then LN₂. If readings are within accuracy specifications noted above, then the temperature circuit on the control is properly calibrated. If not, follow the steps below to calibrate the temperature.

Access the menu by pressing



Press  5 times until the display reads "CALIBRATE".

Press  to enter the calibration mode. The display will read



Submerge the thermocouple into an ice water bath. The number of dots will decrease as the control approaches a good calibration point. 0 to 2 dots will provide a good calibration.

Press  The display will reset.

Submerge the thermocouple into LN₂. The number of dots will decrease as the control approaches a good calibration

point. 0 to 2 dots will provide a good calibration. Press




If the calibration values fall within the expected range, a message "Calibration Successful" will be displayed.

Appendix


External Connector Ratings

Designated Use	Max Rated Voltage/Current Ratings	Connector Type
Power	36 VDC	2.5mm barrel connectors
CAN	N/A	Modular shielded jack
Remote Alarm	300 Volts	5mm terminal block
Solenoid Valve	600 Volts	4.2mm header
Thermocouple	N/A	2 pin thermocouple
Level Sensor	5 amps/Contact	Sealed circular connector

Labels



Main Rating Label



Control & Display Label